


Confectietabel VENOSAN® Prestige 2 en 3 en Prestige 2 Katoen rondbrei beenkousen

maat	small (cm)	medium (cm)	large (cm)	extra large (cm)
	normaal / wijd	normaal / wijd	normaal / wijd	normaal / wijd
cG silb.	48 - 56 / 57 - 66	50 - 60 / 61 - 72	53 - 66 / 66 - 80	57 - 73 / 74 - 90
cG	48 - 66	50 - 72	53 - 80	57 - 90
cF	42 - 58	44 - 62	47 - 68	51 - 77
cE	32 - 41	34 - 44	37 - 48	41 - 54
cD	30 - 37	32 - 39	34 - 42	37 - 47
cC	31 - 39	33 - 41	35 - 44	39 - 48
cB1	24 - 30	26 - 32	29 - 35	33 - 39
cB	19 - 21	21 - 23,5	23 - 27	27 - 32
cY	28 - 34	29 - 35	31 - 37	34 - 40
cA	20 - 24	21 - 25	22 - 27	24 - 30





People.Health.Care.



Confectie assortiment VENOSAN® Prestige 2 en 3 en Prestige 2 Katoen rondbrei beenkousen

	kleur	teenstuk	AB	AD	AD silb. 3 cm/5 cm	AG	AG silb. 5 cm	AGT	AT
VENOSAN® Prestige 2	● beige / ● zwart	open	●	● ●	● ●	●	●	● ●	-
		gesloten	-	● ●	● ●	-	●	-	●
VENOSAN® Prestige 2 Katoen	● beige	open	-	●	●	-	●	●	-
		gesloten	-	●	●	-	-	-	●
VENOSAN® Prestige 3	● beige	open	-	●	●	-	●	●	-
		gesloten	-	-	-	-	-	-	-

AB: enkelkous **AGT:** AG kous met
AD: kniekous heupstuk
AG: dijbeen kous **AT:** panty
silb.: siliconenband

	lengte AD (cm)	lengte AG (cm)	lengte GT voor AT (cm)	voetlengte IA / IZ (cm)
kort	35 - 38			16 - 17 / 22 - 24
normaal	38 - 42	65 - 76	29 - 35	17 - 19 / 23 - 26
lang	42 - 47	76 - 88		19 - 21 / 25 - 28

