




People.Health.Care.

Advieskaart

Confectietabel VENOSAN® Soft 2, Pro 2 en Pro 3 vlakbrei beenkousen

Omvang (cm)	small c-	small	medium c-	medium	large c-	large	x large c-	x large
cD	25 - 30	30 - 35	28 - 33	33 - 38	32 - 38	36 - 42	37 - 43	41 - 47
cC	26 - 31	31 - 36	29 - 34	34 - 39	33 - 39	37 - 43	38 - 44	42 - 48
cB1	21,5 - 24,5	24 - 27	23,5 - 27	26 - 29,5	27 - 30,5	29 - 32,5	29 - 33	31 - 35
cB	19 - 21		21 - 23		23 - 25,5		25,5 - 28,5	
cY	29 - 33		31 - 34		32 - 36		34 - 38	
cA	20,5 - 22,5		21,5 - 23,5		23 - 25		24 - 27	



Lengte (cm)	kort	normaal	lang
ID	36 - 38	39 - 41	42 - 44
IA	16 - 18	17 - 19	18 - 20
IZ	22 - 24	23 - 26	25 - 28



